



# Bariatric+ Patient Chair

[www.knightimaging.co.uk](http://www.knightimaging.co.uk)



**Knight  
Imaging**

EXPERIENCE & EXPERTISE

# Bariatric+ Patient Chair

EXPERIENCE & EXPERTISE

The VELA Bariatric+ Patient Chair provides healthcare staff with a robust solution designed to reduce work-related musculoskeletal disorders (WRMSDs) associated with bariatric patient care including lower back pain, sprained wrists, and shoulder injuries. With an electric lift capability and high weight capacity, it improves the patient's comfort while mitigating the risk of injury to staff.

The Bariatric+ is secure, supportive, and highly manoeuvrable. The chair's easy-running wheels allow staff to move the patient comfortably or allow the patient to position themselves using their legs, encouraging mobility and independence. Supporting a maximum weight of 300 kg, patients can comfortably stand up or sit down using the electric lift function, reducing their reliance on holding onto a nurse or porter, and subsequently reducing the risk of a WRMSD occurring.

Not only does the Bariatric+ reduce the physical strain placed on staff, it is also designed for infection control in wards. It is deliberately designed with smooth surfaces and without recesses to stop the accumulation of dirt. Every part of the chair can be cleaned using common disinfectants including alcohol and chlorine for effective sterilisation.

**PATIENT WEIGHT**  
Up to 300 kgs  
(47 Stone)

Electric lift capability

Easy to clean and disinfect

Easy running wheels

## Specifications

Seat (W x D)	72 x 50 cm
Backrest (W x H)	38 x 43 cm
Seat Height	42 - 57 cm
Armrest Pad (W x L)	9 x 30 cm
Armrest Height (Measured from seat)	0 - 24 cm
Frame Size	58 x 58 cm

# Service and Support

Knight Imaging offers service and support on our range of couches and chairs. We carry out one-off service calls, yearly preventative maintenance visits as well as offering support contracts.

## Planned Preventative Maintenance

Our fully trained engineers service your couches and chairs ensuring that any potential issues are spotted and brought to your attention before they become a problem.

We typically carry out these types of services across entire departments ensuring the best possible lifetime of your couches and chairs. We will always arrange visits in advance to ensure a swift turnaround and to minimise disruption to your ward or department.

A preventative maintenance visit can be arranged whenever you think you need to have your equipment inspected and checked for faults or potential future problems. Just get in touch to discuss your equipment and make a booking.

## Breakdown Response

Our experienced engineers are available to diagnose and repair in times of emergency.



**Knight Imaging Limited**  
Innovation House, 6 Europa View  
Sheffield, S9 1XH, UK  
**E-mail:** [sales@knightimaging.co.uk](mailto:sales@knightimaging.co.uk)  
**Tel:** +44 (0)114 267 0482  
[www.knightimaging.co.uk](http://www.knightimaging.co.uk)

